



EVO

 Visian ICL™



Evolution in Visual Freedom.™



The EVO Visian ICL Advantages

Many vision correction procedures promise an improved level of vision, but few vision correction alternatives offer the quality and features found with the EVO Visian ICL.

These include:

- Sharp, Clear Vision^{12,3}
- Does not induce dry eye syndrome⁴
- Excellent Night Vision⁵
- Great for those with thin corneas^{6,7}
- UV protection built into the biocompatible lenses
- A typically quick, 20-30 minute procedure and most people are able to resume daily activities in just a few short days with clearer vision
- Treats moderate to high myopia (nearsightedness), Hyperopia* (farsightedness) and astigmatism.

* Hyperopic lenses require an additional tiny opening in the colored portion of your eye to ensure that the fluid inside your eye flows properly. Consult your physician for more details.



What It Is and How It Works

The EVO Visian ICL is an Implantable Collamer Lens that works with the eye to correct vision. Unlike traditional contact lenses that go on the surface of the eye, the EVO Visian ICL is positioned inside the eye between the iris (the colored part of the eye) and the natural lens where it stays indefinitely. If your vision changes dramatically, however, it can be removed from the eye.



Technologically Advanced Lens Material

Unlike other lenses, the EVO Visian ICL is made from Collamer, a technologically advanced lens material exclusive to STAAR Surgical. Because of the collagen in Collamer, it is biocompatible with the eye (not seen as a foreign object to the body) and provides exceptional quality of vision.¹ The lens rests in position while accurately correcting vision.¹



Designed to be Permanent, Yet is Removable

One of the most common concerns people have about laser based correction is that the procedure permanently changes the shape of their eye by removing corneal tissue. An advantage of the EVO Visian ICL is that it can permanently correct your vision, yet no natural tissue is removed. If for the unlikely event you need or wish to remove your EVO Visian ICLs in the future, your doctor can accommodate you.



Harmony with your Natural Eyes

For people who value health and wellness, the EVO Visian ICL could be the choice for you. The Collamer lens material is biocompatible, meaning it is designed to function in harmony with your natural eyes.



Legacy of Success

The EVO Visian ICL has been approved in Europe since 1997 and in the U.S. since 2005. In that time, over 550,000 eyes have been implanted with Visian ICLs by certified, specialized surgeons.⁸ In a clinical trial, patients reported a 99% satisfaction rate.⁹



Treatment Flexibility for the Future

Your vision changes as you grow older and eventually almost everyone will need cataract surgery. Because EVO Visian ICL preserves the structure of the eye, the lens can be removed and replaced if your vision changes substantially, or if you have another procedure at any time.



Premium Outcomes¹⁰ with Exceptional Patient Satisfaction⁹ and Established Safety¹¹

Why do people who demand premium vision performance select the EVO Visian ICL? Because they choose a technology to meet their individual vision needs. There have been over 550,000 patient eyes treated with the Visian ICL worldwide⁸ with a satisfaction rate of 99%.⁹



Established Safety Profile

A study has shown the EVO Visian ICL has an established safety profile. Studies have shown that the EVO Visian ICL provides safe postoperative results with very high patient satisfaction of over 99%.⁹





Good For You to Know

- The Vision ICL is approved for use in adults from 21 to 45 years of age:
- With the correction of myopia in adults ranging from -0.5 D to -18.0 D with or without astigmatism up to 6.0 D
- With the correction of hyperopia in adults with hyperopia ranging from +0.5 D to +10.0 D with or without astigmatism up to 6.0 D



Results may vary and the EVO Visian ICL may not be right for everyone. Like all surgical treatments of the eye the EVO Visian ICL carries risks including: need for a second surgery, temporary swelling or inflammation, cataracts, increases in eye pressure, retinal detachment, endothelial cell loss, and reduced visual acuity. Consult with your healthcare provider to learn more about the EVO Visian ICL, possible complications, and to see if you will be a candidate for the EVO Visian ICL.

The EVO Visian ICL should not be used in people who have the following characteristics

- Anterior Chamber Depth (ACD) of $< 3.0\text{mm}$
- Anterior chamber angle $<$ Grade III as determined by gonioscopic examination
- Are pregnant or nursing
- Patients with low/abnormal corneal endothelial cell density, Fuch's dystrophy, or other corneal pathology
- Ocular hypertension in either eye.
- Any cataract in the operative eye or nontraumatic cataract in the fellow eye.
- Previous or pre-existing ocular disease that would preclude post-operative visual acuity of 20/60 or better.
- Patients who are amblyopic or blind in the fellow eye.
- Keratoconus
- Primary open angle or narrow angle glaucoma

A person with blonde hair, wearing a green backpack and teal shorts, is seen from behind, walking through a cave. They are holding a trekking pole. A large, weathered log is positioned diagonally across the cave opening, acting as a bridge. The cave walls are made of layered, reddish-brown rock. The floor is wet and reflective. The scene is dimly lit, with light coming from the cave entrance.

Time to move on from contacts and eyeglasses

Are contacts and eyeglasses a constant hassle and a huge expense?

The EVO Visian ICL can give you the freedom to live a life uninhibited from the day-to-day lifestyle hassles of glasses and contacts. Your favorite activities and new adventures are on the horizon in usually a few short days!